



GYNIUS 2 SALI QUATERNARI

SANITIZING CONCENTRADED DETERGENT

FEATURES AND USE

It is a detergent assembled base sanitizing of quaternary salt of ammonium, agent sanitizing to fit ample spectrum of action on varied hard surfaces. Particularly suitable for the daily maintenance and the rapid sanitation of rolled surfaces, you lacquer, in ceramics and floors in gres, marble granulated etc., in environments what hospitals, nursing homes, cafeterias and collectivity. Pleasant scent of mentholated pine.



METHOD OF USE & DILUTION:

After opportune dilution through the system of dilution GYNIUS, is usable both manually with the employment of the mop that with cars washer and dryer. After the dilution and the use of the gotten detergent, to rinse e/o to always dry the treated surfaces.

PACK SIZE:

Jerrycan 10/20Kg - Unit-doses 15ml – 160 pcs/carton - Unit-doses 30ml – 120 pcs/carton

TECHNICAL DATA:

Physical state: liquid

Colour: bright blue

Odour: Balsamic pine

pH: 8,66 +/- 0,50

Specific weight: 1,020 +/- 0,006 g/cm³

COMPOSITION:

See safety data sheet. Ingredients compliant with CE Reg. 648/2004

CAUTIONS:

Causes severe skin burns and eye damage. Very toxic to aquatic life.

May produce an allergic reaction.

Precautionary statements:

Avoid release to the environment. Wear protective gloves / protective clothing / eye protection / face protection.

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON CENTER or doctor / physician.

Contains: Quaternary ammonium compounds, Benzyl-C12-16-alkyldimethyl, chlorides BRANCHED ALKYL ALCOHOL (C11-C13) ETHOXYLATES.



Warning



AR-CO CHIMICA s.r.l. Via Canalazzo, 22/24 – 41036 MEDOLLA - ITALY
Tel 0535.58890 – info@arcochimica.it – www.arcochimica.it

PROFESSIONAL USE