




GYNIUS 4 CITRICO

CONCENTRATED CITRIC ACID-BASED DETERGENT
FOR SANITARY FITTINGS



FEATURES AND USE

GYNIUS 4 WITH PEACH/FLOWER SCENT, is a citric acid-based maintenance detergent suitable for daily cleaning of washbasins, taps, sanitary fittings, showers and ceramic tiles. Removes calcareous residues mixed with soap. Fast drying product with good glossing characteristics that makes daily cleaning operation easy, and leaves a pleasant and long-lasting scent. Used following the proper method of use it is delicate on metal surfaces.



METHOD OF USE & DILUTION:

Dilute GYNIUS 4 with GYNIUS dilution system and spray the solution on the surface to be cleaned; rinse and/or dry always the treated surfaces with cloth or paper.

METHOD OF USE GYNIUS PACK : Pour the 100 ml unit dose product into the special Gynius Pack bottle, add 500 ml of water to obtain 600 ml of cleaning solution ready to use.

PACK SIZE:

Jerrycan 10/20Kg – Unit-doses 100ml – 60 pcs/carton

TECHNICAL DATA:

Physical state : Liquid

Colour: Brilliant Pink

Scent: Peach/Flower

pH as is 2,8+/-0,5

Specific weight: 1,024+/-0.009 g/cm³

COMPOSITION:

See safety data sheet. Ingredients compliant with CE Reg. 648/2004

CAUTIONS:

Causes serious eye damage. Causes skin irritation.

Precautionary statements:

Wear protective gloves / protective clothing / eye protection / face protection.

IF ON SKIN: Wash with plenty of soap and water. Immediately call a POISON CENTER or doctor/ physician.

Contains: ALCOHOL ETHOXYLATED C12/15 >6-9EO Methanesulfonic Acid.



Warning



AR-CO CHIMICA s.r.l. Via Canalazzo, 22/24 – 41036 MEDOLLA - ITALY
Tel 0535.58890 – info@arcochimica.it – www.arcochimica.it

PROFESSIONAL USE