



# GYNIUS CANDOR

POWDER DETERGENT WITH PERBORATE FOR  
LAUNDRY WASH



## FEATURES AND USE

Gynius CANDOR is an atomized detergent in unit-dose sachets for laundry machine wash. Gynius CANDOR can be used with any type of laundry, domestic and industrial. The high content of sequestrants allows using Gynius CANDOR with any degree of water hardness, in case of very hard water the dose must be increased. Gynius CANDOR does not produce calcareous residue on laundry or on the machine, thus helping to maintain the fabric soft. Gynius CANDOR contains anionic, non-ionic surfactants, phosphates, phosphonates, sodium soap, stabilizers, anti-foaming agent, colorants and fragrance.

### METHOD OF USE & DILUTION:

Wash in washing machine : put the sachet of Gynius Candor into the washing compartment.

#### PREWASHING:

SOFT WATER ( 0-15° F): 8 gr for every kg of dry fabric

SLIGHTLY HARD WATER (15-25°F): 15 gr for every kg of dry fabric

HARD WATER (more than 25°F): increase dosage of 1 gr for every hardness grade above 25° F.

#### WASHING:

SOFT WATER ( 0-15° F): 18 gr for every kg of dry fabric

SLIGHTLY HARD WATER (15-25°F): 25 gr for every kg of dry fabric

HARD WATER (more than 25°F): increase dosage of 1 gr for every hardness grade above 25° F.

### PACK SIZE:

Unit-doses 60g – 30 pcs/carton

### TECHNICAL DATA:

Physical state : powder

Colour: White

Odour: TYPICALLY

pH: 10,5 -11

Specific weight :720kg/lt+/- 0,05kg/lt

### COMPOSITION:

See safety data sheet. Ingredients compliant with CE Reg. 648/2004

### CAUTIONS:



Warning



AR-CO CHIMICA s.r.l. Via Canalazzo, 22/24 – 41036 MEDOLLA - ITALY  
Tel 0535.58890 – info@arcochimica.it – www.arcochimica.it

**PROFESSIONAL USE**